



Welch Allyn

ProBP™ 2000 Digital Blood Pressure Device




Fast and accurate automated blood pressure readings at your fingertips



WelchAllyn®

Welch Allyn ProBP 2000—the IDEAL BLOOD PRESSURE SOLUTION for your care setting

We've combined two of our innovative technologies to enhance exams for you and your patients.



Welch Allyn SureBP® technology

Our SureBP inflation-based algorithm is used across Welch Allyn vital signs devices, helping clinicians in every care setting capture automated, fast and accurate blood pressure readings.

- Blood pressure readings are captured more comfortably on inflation and take about 20 seconds
- Meets AAMI standards for blood pressure accuracy
- One-button touch triggers the device to start capturing on cuff inflation
- 97% completion rate in the presence of motion¹

“SureBP technology offers great advantages in speed and patient comfort... clinicians can have added confidence in this new technology².”

Dr. Bruce S. Alpert | University of Tennessee,
AAMI Sphygmomanometer Committee Chair

Welch Allyn FlexiPort® Blood Pressure Cuff technology

Using the correct cuff size is an important factor in obtaining accurate readings. Welch Allyn FlexiPort cuffs come in a wide variety of sizes to help address a larger range of your patient population, helping you ensure the right cuff size is used each time, to support accuracy.

- Unique “snap-fit”, single-point connection helps remove the hassle of traditional cuff tubes and connectors, to better support a busy workflow
- The ProBP 2000 is compatible with FlexiPort Reusable cuffs, ranging from child up to large adult and thigh (15-55 cm range)
- When combined with SureBP technology, helps capture readings on inflation in about 20 seconds¹



When looking for a clinical digital blood pressure device, **accuracy** and **reliability** may be your biggest concerns. Welch Allyn has designed the new ProBP 2000 digital blood pressure device with technology features for results you can feel confident in.



ProBP 2000 Digital Blood Pressure Device Product Specifications

External Dimensions	Width: 94 mm (3.70 in) Length: 142 mm (5.59 in) Height: 66 mm (2.59 in)
Weight	283 g (0.62 lbs) without batteries
Display	Digital LCD V.A. 68 mm x 90 mm, pressure and pulse displayed
Power Options	4 AA alkaline batteries, 6VDC (included) AC power supply (107041) Input: 100-240 V, 50-60 Hz, 400 mA Output: 6 V, 1
Method	Oscillometric
Measurement Range	Rated cuff pressure: 0 mmHg - 299 mmHg (0 kPa - 40 kPa) Measurement pressure: Systolic: 50-260 mmHG Diastolic: 25-220 mmHG Pulse: 40 to 199 beats/minute
Accuracy	Pressure: +/- 0.4kPa (3 mmHg) Pulse Value: +/- 4%

Parts & Accessories*

CONFIGURATIONS

- 2000** ProBP 2000 Digital Blood Pressure Device, size 11 Adult FlexiPort Reusable BP cuff, 4 AA batteries
2000-P ProBP 2000 Digital Blood Pressure Device, size 11 Adult FlexiPort Reusable BP cuff, 4 AA batteries and power supply

ACCESSORY

- 107041** Power Supply
(includes US, EU, UK,
AU power adapters)

FLEXIPORT BLOOD PRESSURE CUFFS

- | | |
|--|--|
| REUSE-09 Child (15-21 cm) | REUSE-12 Large Adult (32-43 cm) |
| REUSE-10 Small Adult (20-26 cm) | REUSE-12L Large Adult Long (32-43 cm) |
| REUSE-11 Adult (25-34 cm) | REUSE-13 Thigh (40-55 cm) |
| REUSE-11L Adult Long (25-34 cm) | |

*Check with your local Welch Allyn representative for a full list of configurations and accessories available in your region.

To learn more, visit emeai.welchallyn.com or contact your Welch Allyn sales representative.

¹ Welch Allyn Motion Tolerance Clinical Study, August 2016

² Alpert, Bruce S., Clinical evaluation of the Welch Allyn SureBP® algorithm for automated blood pressure measurement. Blood Pressure Monitoring, Vol. 12 No. 4, August 2007.



Welch Allyn EMEA Customer Care
IDA Business Park - Johnstown
Navan, County Meath, Ireland
(p) +353 46 906 7700
europe@welchallyn.com

EMEAI.WELCHALLYN.COM

WelchAllyn®

Hill-Rom®

©2018 WELCH ALLYN, INC. ALL RIGHTS RESERVED. DOC. NO: SM4185EU REV.A MC14918, MAY 2018

The ProBP 2000 Digital Blood Pressure Device is intended for use in measuring blood pressure and heart rate in pediatric and adult patient populations. Medical device (93/42/EEC); Class IIa. Notified Body: TÜV Rheinland LGA Products GmbH (CE0197) Manufacturer: Guangdong Transtek Medical Electronics Co., Ltd. No. 105 Dongli Road, Torch Development District, Zhongshan, 528437, Guangdong, China. EC Rep.: MDSS-Medical Device Safety Service GmbH, Schiffgraben 41, 30175 Hannover, Germany. This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Welch Allyn recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices. Welch Allyn reserves the right to make changes without notice in design, specifications and models. The only warranty Welch Allyn makes is the express written warranty extended on the sale its products.