



4302 Diaper and Organ Scale

Durable. Long Lasting. Designed for Demanding Healthcare Environments.

- With increments of one gram, you can measure urine output
- Increments of grams provides one-to-one relation to cc's of urine for easy, rapid calculations
(Pounds and ounces may also be selected.)

Configurations

4302-XX-X	Standard capacity - 3,000g (X), standard weight (X) and battery power (X)
4302-XK-X	Standard capacity - 3,000g (X), grams only (K) and battery power (X)
4302-XX-B	Standard capacity - 3,000g (X), standard weight (X) and US AC power (B)
4302-XK-B	Standard capacity - 3,000g (X), grams only (K) and US AC power (B)
4302-EX-X	Extended capacity - 5,000 g (E), standard weight (X) and battery power (X)
4302-EK-X	Extended capacity - 5,000 g (E), grams only (K) and battery power (X)
4302-EX-B	Extended capacity - 5,000 g (E), standard weight (X) and US AC power (B)
4302-EK-B	Extended capacity - 5,000 g (E), grams only (K) and US AC power (B)

Specifications

Weighing Capacity	Standard capacity: 3,000 g Extended capacity: 5,000 g
Dimensions	8" w x 11" d x 3" h
Platform	8.5" w x 8.5" d
Resolution	1 g
Readout	Digital LED in grams, pounds or ounces. Automatic shutoff.
Power Source	6 C-sized alkaline batteries. Power adapter available.
Weight	9 lbs
Warranty	One year
Optional Accessories	Power adapter available
Automatic Zero	N/A
Weight Recall	N/A
Reweight	N/A
Computer Capability	N/A

Contact us at scale-tronix@welchallyn.com or visit www.welchallyn.com to learn more.

Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3694