


Document Description: TSB –VS100CR - IR Filter	Document Number: 20014255 Version: A
	
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only.	

<h1>Technical Service Bulletin</h1>	
Product: SPOT VISION SCREENER	Date: 2016-09-23
Subject: IR Filter cost reduction replacement	
HW Version(s) Affected: VS100CR	SW Version(s) Affected: VSD 3.0.04.06 and later
Serial Numbers Affected: N/A. Running change	Lot or Date Code Affected: N/A. Running change
Classification: Informational Only	
Distribution: <input type="checkbox"/> Customer Care <input checked="" type="checkbox"/> Product Service <input type="checkbox"/> Field Service <input type="checkbox"/> ASPs <input type="checkbox"/> Distributors <input type="checkbox"/> Customers <input type="checkbox"/> Company Confidential	
<u>THIS DOCUMENT CONTAINS CONFIDENTIAL INFORMATION</u>	
It can be distributed or made available to only those individuals, companies, and organizations which have current Welch Allyn Confidentiality Agreements, and have a need to know.	
Training Required: <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	
Summary: As part of the DMACs cost reduction effort, IR filter on the SPOT Vision Screener is replaced by a cost reduced equivalent. This is a running change per ECR 400132093. Current IR Filter p/n 20020 will be replaced by p/n 724385 on 30028_RBOM. The following service procedures were updated to remove reference to the specific part numbers and thus leaving only part description: <ol style="list-style-type: none"> 1. DIR #20013679 version D updated to version E – on the 30028_RBOM 2. DIR #20011080 version D updated to version E – on the VS100_RBOM 	

** Training is required if the document includes information on how to change/repair/fix the product - and only applies to TSBs

Version	Sec, Pg, Para Changed	Change Made	Date Version Created	Version Created By (initials)
A	N/A	Initial Release	2016-09-23	SC