

## Welch Allyn Video Otoscope Comparison - Human Version

	<b>CVO 23100</b>	<b>Digital MacroView 23920</b>
Signal Output	Analog	Digital
Video Image Display Requirements	RCA/Composite, TV/Video monitor 9 - 13"	Computer w/Windows XP Service Pack 2.0 or higher
Image Resolution	510 x 492	1280 x 1024
Power Sources Supported	3.5v Rechargeable handles only	3.5v Rechargeable handles & Wall Transformer
Power to Device	Handle supplies camera and lamp	Handle supplies lamp, computer supplies camera
Pneumatic Port	No	Yes
Focus	No - Fixed Focus	Yes - Variable
Image Capture	Not as sold. Requires frame grab software and USB video bus. Images are degraded from transfer from analog to digital	Yes - megapixel still images
Zoom	No	Yes
Image Flip	No	Yes
Video Clip Tool	No	Yes
08 List Price Head only	\$1,265.45	\$1,185.00
08 List Price w/ 71000-A handle	\$1,220.35	\$1,346.00