

| <b>Model</b> | <b>Description</b>      | <b>Diopters</b>                  | <b>Focal Length</b> |
|--------------|-------------------------|----------------------------------|---------------------|
| 20000/25020  | 3.5V Otoscope           | 8.8                              | 114mm               |
| 20100/20200  | Standard Otoscope       | 10                               | 100mm               |
| 20160/20260  | Vet Otoscope            | 8                                | 125mm               |
| 20151/20251  | 12 Diopter Otoscope     | 12                               | 83mm                |
| 21110        | 2.5V Pocketscope        | 9.8                              | 102mm               |
| 21600/21700  | Operating Otoscope      | 7.9                              | 127mm               |
| 21660/21760  | Vet Operating Otoscope  | 6.5                              | 152mm               |
| 20202        | Consulting Otoscope     | 6.5                              | 153mm               |
| 20262        | Vet Consulting Otoscope | 5.6                              | 177mm               |
| 23540        | Mag Lens Otoscope       | 2                                | 500mm               |
| 24031        | 2.5V Otoscope           | 8.8                              | 114mm               |
| 29090        | Otoscope                | 9.8                              | 102mm               |
| 23810, 23820 | MacroView Otoscope      | Diopter adjustment is continuous |                     |

| Magnification  |
|----------------|
| 2.2X           |
| 2.5X           |
| 2.0X           |
| 3.0X           |
| 2.5X           |
| 2.0X           |
| 1.7X           |
| 1.7X           |
| 1.4X           |
| 0.5X           |
| 2.2X           |
| 2.5X           |
| 3.6X at 27.6mm |