





















# FLEXIPOINT™ TUBE SET REFERENCE CHART

	Part Number	Fitting(s)/Tubing/Adapter	Unit of Measure	Tube	Connector(s)	Monitor Compatibility
<b>One-Tube</b>	1-BARB		10	(1) 8"	Double hose barb (5082-171)	This is used to semipermanently connect the FlexiPort fitting to any one-tube device. It requires removing the original fitting from the device tube.
	1-SC		10	(1) 8"	Screw (5082-164)	One-tube screw-type devices (e.g., Welch Allyn Spot, Welch Allyn VSM, Welch Allyn Atlas)
	1-MQ		10	(1) 8"	Locking (5082-191)	Welch Allyn one-tube locking-type devices
	1-HP		10	(1) 8"	Bayonet (5082-190)	One-tube Bayonet-type devices (e.g., HP, Philips, Spacelabs, DataScope)
	1-TPS		10	(1) 24"	Tri-purpose (5082-168)	One-tube hand aneroid devices
	1-TPSE		10	(1) 24"	Tri-purpose (5082-168) Bagged	European hand aneroid devices
	1-TPL		10	(1) 31"	Tri-purpose (5082-168)	One-tube hand aneroid devices
1-TPLE	10		(1) 31"	Tri-purpose (5082-168) Bagged	European hand aneroid devices	
<b>Two-Tube</b>	2-BARB		10	(2) 8"	Double hose barb (5082-171)	This is used to semipermanently connect the FlexiPort fitting to any two-tube device. It requires removing the original fittings from the device tubes.
	2-SC		10	(2) 8"	Screw (5082-164)	Two-tube screw-type devices (e.g., Critikon Dinamap)
	2-SCC		10	(2) 8"	Screw (5082-164) Clamp on tube	Two-tube screw-type devices (e.g., Critikon Dinamap). When Clamp is activated, this can be used with one-tube devices (e.g., Welch Allyn Spot, Welch Allyn VSM, Welch Allyn Atlas)
	2-MQ		10	(2) 8"	Locking (5082-191)	Two-tube locking-type devices (e.g., Spot LXi, GE Marquette)
	2-MF		10	(2) 8"	Locking (Male 5082-172 and Female 5082-191)	Two-tube mated locking-type devices (e.g., Datex Ohmeda)
	2-TPLE		10	(2) 31"	Tri-purpose (5082-168) Bagged	European wall, mobile, and pocket aneroids
	2-TPSE		10	(2) 24"	Tri-purpose (5082-168) Bagged	European wall, mobile, and pocket aneroids
	2-TP		10	(1) 8" (1) 13"	Tri-purpose (5082-168) Bagged	Wall, mobile, and pocket aneroids (e.g., 767, 509, DS48A)
	2-SRC		10	(1) 8" (1) 13"	Tri-purpose (5082-168)	Wall, mobile, and pocket aneroids (e.g., 767, 509, DS48A)

	Part Number	Fitting(s)/Tubing/Adapter	Unit of Measure	Tube	Connector(s)	Monitor Compatibility
<b>Two Tube with Bulb</b>	2-BVL		10	(1) 8" (1) 13"	Tri-purpose (5082-168) Large Bulb	Wall, mobile, and pocket aneroids (e.g., 767, 509, DS48A).
	2-BVLC		10	(1) 8" (1) 13"	Tri-purpose (5082-168) Large	Wall, mobile, and pocket aneroids. (e.g., 767, 509, DS48A) When clamp is activated it can work with one-tube aneroids (e.g., DS58, DS66)
	2-BVS		10	(1) 8" (1) 13"	Tri-purpose (5082-168)	Wall, mobile, and pocket aneroids (e.g., 767, 509, DS48A)
	2-BVSC	10	(1) 8" (1) 13"	Tri-purpose (5082-168)	Wall, mobile, and pocket aneroids (e.g., 767, 509, DS48A). When clamp is activated it can work with one-tube aneroids (e.g., DS58, DS66)	
	2-BVD		10	(1) 8" (1) 13"	Tri-purpose (5082-168)	Wall, mobile, and pocket aneroids (e.g., 767, 509, DS48A).
	2-BVDC		10	(1) 8" (1) 13"	Tri-purpose (5082-168) Econo Bulb & Clamp"	Wall, mobile, and pocket aneroids (e.g., 767, 509, DS48A). When clamp is activated it can work with one-tube aneroids (e.g., DS58, DS66)
<b>Other Accessories</b>	5082-188		10	8" FlexiPort tubing		
	5082-189		10	13" FlexiPort tubing		
	PORT-1		10	One-tube FlexiPort fitting without any tubes		
	PORT-2		10	Two-tube FlexiPort fitting without any tubes		
	5082-159		1	Wash cap that teathers to FlexiPort on cuffs		This allows users to close off the FlexiPort during washing. Also can be used to prevent debris from entering the port.
	1SRCL		1	DuraShock Classic Adapter with one 24" tube and a large bulb with premium valve		For DS48 Integrated Gauges to allow them to work with FlexiPort cuffs.
	1SRCS		1	DuraShock Classic Adapter with one 24" tube and a small bulb with premium valve		For DS48 Integrated Gauges to allow them to work with FlexiPort cuffs.
	1DSAL		1	DuraShock Adapter with one 24" tube and a large bulb with premium valve		For DS45 Integrated Gauges to allow them to work with FlexiPort cuffs.
	1DSAS		1	DuraShock Adapter with one 24" tube and a small bulb with premium valve		For DS45 Integrated Gauges to allow them to work with FlexiPort cuffs.
	1DSLCL		1	DuraShock Adapter with one 24" tube and a economy bulb & valve"		For DS44 Integrated Gauges to allow them to work with FlexiPort cuffs.