

IHE Integration Statement

| Vendor | Product Name | Version | Date |
|-------------|--------------|---------|-------------|
| Welch Allyn | Connex® VM | v2.0 | 20-Jan-2012 |

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

| Profiles Implemented | Actors Implemented | Options Implemented |
|--|-------------------------------|---|
| Consistent Time | Time Client | Secure NTP. |
| Device Enterprise Communication (2011) | Device Observation Reporter | |
| Patient Administration Management | Patient Demographics Consumer | Merge |
| | Patient Encounter Consumer | Inpatient / Outpatient Encounter Management |

| | | |
|--|---|--|
| Internet address for vendor's IHE information: http://www.welchallyn.com/documents/Patient%20Monitoring/Vital%20Signs%20Capture/Connex%20Data%20Management%20System/IHE%20Integration%20Statement.pdf | | |
| Links to Standards Conformance Statements for the Implementation | | |
| HL7 | http://www.welchallyn.com/products/en-us/x-18-mu-96-1179856474703 . | |
| DICOM | | |
| Links to general information on IHE | | |
| In North America: www.ihe.net | In Europe: www.ihe-europe.net | In Japan: www.ihe-j.org |